The Liquid Jet Characterization Pipeline at SLAC: Advancing X-Ray Laser Experiments

Dimitra Manatou (1,2) – dmanatou@slac.stanford.edu, B. Sublett (1), Y. Willard (1), W. P. Chang (1), C. Hampton (1), K.Nicolay (1), K.Karpos (2), R.Kirian (2)

1. SLAC National Accelerator Laboratory, USA 2. Department of Physics, Arizona State University, USA







INTRODUCTION

The liquid jet characterization pipeline, originally developed at Arizona State University (ASU), was replicated at SLAC National Accelerator Laboratory during this internship. This adaptation helps in evolving the experimental tools available for sample delivery injector characterization in crystallography and solution scattering experiments. The transition not only entails the physical duplication of the setup but also incorporates the use of proven Python-based software to ensure methodological consistency and expand experimental possibilities.

At SLAC, enhancements like the integration of LEDs for improved visualization and measurement extend the utility of this pipeline. It can serve as an asset for LCLS staff, enabling them to efficiently test injectors and conduct miniexperiments. The automated data collection and image processing capabilities of this system accelerate the production of diagrams that correlate flow conditions with jet properties, optimizing jet configurations to meet precise experimental demands. This evolution not only boosts experimental efficiency but also reduces the risk of

SOFTWARE






Figure 3: The Odysseus interface for controlling GDVNs and collecting stacks of jet images. Bronkhorst gas control and HPLC liquid flow control Python examples are seen below.

EXPECTED RESULTS





expensive failures during pivotal research phases.

SCHEMATICS





Figure 1: (A) Dual pulse triggering scheme. LED flashes are 1us apart.



EXPERIMENTAL SETUP




Figure 4: LED lighting up the vacuum chamber with the sample delivery injector in order to image it through the software shown above.



Figure 6 (above): Phase diagram indicating jet speed for each "control knob" setting (liquid volumetric flow rate and gas mass flow rate). Marker size indicates jet diameter.

FUTURE STEPS

Data Collection at the ICL Testing Station: The software has been adapted to efficiently collect and analyze data at the Injector Characterization Lab (ICL). Implementation in ASC Wet Labs: To support cryo-electron microscopy (cryoEM) experiments, the software was integrated into the ASC wet labs. Potential Integration into Experimental 3 Hutches: Plans are underway to integrate the software into one of the experimental hutches. This integration is envisioned as a method for continuous monitoring and control of the sample delivery systems, ensuring optimal operation during critical experimental runs.




Figure 7: ASC set up





Figure 2 : Schematic of the triggering and set up of electronics and pulses.






Figure 7: CXI hutch.

ACKNOWLEDGEMENTS

Use of the Linac Coherent Light Source (LCLS), SLAC National Accelerator Laboratory, is supported by the U.S. Department of Energy, Office of Science, Office of Basic Energy Sciences under Contract No. DE-AC02-76SF00515.

Figure 5: ICL set up

MEC end-station and itsparticipation in LaserNetUS is supported by the U.S. Department of Energy, Office of Science, Office of Fusion Energy Sciences under Contract No. DE-AC02-76SF00515.

REFERENCES

- 1. 3D printing of gas-dynamic virtual nozzles and optical characterization of high-speed microjets, Nazari, R., et al, 2020
- 2. Femtosecond X-ray protein nanocrystallography, Chapman, H., et al, 2011
- 3. Generation of Steady Liquid Microthreads and Micron-Sized Monodisperse Sprays in Gas Streams, Calvo, G., 1998
- 4. Gas Dynamic Virtual Nozzle for Generation of Microscopic Droplet Streams, DePonte, D., et al, 2008
- 5. High-resolution protein structure determination by serial femtosecond crystallography, Boutet, S., et al, 2012