


CLEAR FOR LAUNCH:

Standardizing Cellulose and Pluronic Protocols for Crystal Delivery

Ava Markert^{1,2}, Henry Li¹, Kathryn Olsen¹, Raymond G. Sierra¹

1. LCLS, SLAC National Accelerator Laboratory, USA. 2. Rensselaer Polytechnic Institute, Troy, NY, USA.



BACKGROUND

XFEL crystallography needs a steady stream of protein crystals delivered without clumping or settling

- Crystals suspended in a viscous “carrier” gel
- Extruded through a narrow capillary.
- Increased viscosity → continuous liquid → reduces sample consumption.

WHY CELLULOSE AND PLURONIC

- Why hydrogels?



Figure 1.1. Main characteristics of hydrogel use in XFEL.

- Hydroxyethyl cellulose (HEC) and Pluronic F-127 are hydrophilic carrier media, alongside LCP.
- Both gels are much simpler to prepare than LCP.

METHODS

HEC gels:

- Concentrations tested: 10–18% w/v
- Tested with and without heat during mixing
- Standardized mixing procedure for consistency




Figure 2.1 Flowchart diagram of hydroxyethyl cellulose protocol.

Pluronic F-127:

- Screened 20–40% w/v (5% steps)
- 20% did not solidify and 40% did not dissolve
- Narrowed to 25–30% w/v (1% steps)




Figure 2.2 Flowchart diagram of Pluronic F-127 protocol.

Both gels:

- ~15 min active mixing
- Minimum 2-day rest before use
- Extrusion tested via HPLC pump + high viscosity capillary extruder

RESULTS

HEC

- No heat → solid layer at the bottom, loose at the top (poor homogeneity)
- With heat → more homogeneous, but higher concentrations solidified faster.
- Heating is more homogeneous, and it is recommended to use 12% for future user experiments.




Figure 3.1 Image of both heated and nonheated HEC gels ranging from 10–18% (left to right) concentration on a 96-well plate. Rows A-D are heated HEC, and rows E-H are non-heated.

Pluronic F-127

- Faster, more complete dissolution with revised protocol




Figure 3.2 Image of Pluronic F-127 gels ranging from 25–30% concentration.

Extrusion

- Both gels extruded successfully through the high-viscosity extruder.
- There was a slight problem with the gel strands curling while they were being extruded at lower flow rates.




Figure 3.3. Diagram of a high-viscosity extruder.

DISCUSSION

- Recommendation: 11–13% HEC → most homogeneous, fewest trapped bubbles, and viscous.
- Pluronic was homogeneous between 25–30% but looser compared to the cellulose.
- Limitations: evaluation was mostly qualitative

Next steps:

- Quantify viscosity via viscometer
- Test in XFEL

CONCLUSION

- Standardized HEC and Pluronic F-127 protocols improved observed batch consistency and homogeneity
- Recommended 11–13% HEC with heat and 25–30% Pluronic as the starting condition for future users.
- Both gels are far simpler to prepare than LCP.

ACKNOWLEDGEMENTS

I would like to thank my mentor, Ray Sierra, for his guidance and support throughout this project. I am also grateful to Kathryn Olsen for her help with previous procedures, Brigid Baker for helping me get adjusted at SLAC, and Henry Li for his assistance with the methods process. Finally, thank you to SLAC National Accelerator Laboratory and LCLS for providing the facilities and opportunity to carry out this work.

REFERENCES

[1] Schmolka, I.R. Artificial skin I. Preparation and properties of Pluronic® F-127 gels for treatment of burns. *J. Biomed. Mater. Res.* 6, 571–582 (1972). <https://doi.org/10.1002/jbm.820060609>
[2] Kováčsová, G., Grünbein, M.L., Kloos, M. *et al.* Viscous hydrophilic injection matrices for serial crystallography. *IUCrJ* 4, 400–410 (2017). <https://doi.org/10.1107/S2052252517005140>
[3] Benyounes, K., Remli, S. & Benmounah, A. Rheological behavior of Hydroxyethylcellulose (HEC) solutions. *J. Phys.: Conf. Ser.* 1045, 012008 (2018). <https://doi.org/10.1088/1742-6596/1045/1/012008>